

STELLANT

CT Injection Systems

For four decades, Medrad has been changing the way medicine is practiced, by developing some of the most technologically advanced medical devices in healthcare. Now we continue our tradition of contrast delivery innovation with the Stellant Sx Injection System. Stellant Sx is the result of listening to and meeting the needs of healthcare professionals in the most demanding CT environments.

Medrad's Stellant Sx is the Right Choice for the CT Environment of Today and Tomorrow

/Prime Button

• Pushes a predetermined amount of contrast into the tubing to displace air.

Remote Check For Air

• Confirm that no air is in fluid path at the injector head to reduce risk of an air injection.

/Remote Arming

meorao

• Arm the injection system at the injector head.



Single Xpanded

Improved Syringe Design Saves Time and Increases Throughput Capability

Quick and Easy Syringe Installation

- Simply snap the syringes on in any orientation.
- Plunger advances automatically if contrast syringe is empty.
- Just twist to remove the syringe following the procedure plunger will retract automatically into position for the next case.

Integral Auto Load

• Simply connect the fluid path, and the syringe is automatically loaded with the volume of contrast required for the protocol.



Priming Tube

Mini reservoir captures contrast or saline overflow while priming patient tubing.

The New Stellant User Interface Provides Many Powerful Features

The Stellant Interface was designed from the ground up to be the most intuitive and user-friendly system available.



Color Touch Screen

- State-of-the-art color touch screen makes all system operations more intuitive and easier to view.
- Easy-to-read, anti-glare screen.
- Move easily from page-topage without scrolling.
- Stores and recalls up to 32 different protocols.



was



Pressure Monitor Graph

- Displays the pressure being developed in the syringe during an injection.
- Allows operator to reduce patient risk by providing real-time visual graph of pressure.

Flow Profile

• Provides a real-time graphic display of the entire injection profile. You will always know at a glance what has been programmed, and where the injection system is in the delivery of the protocol.

Stop Watch Feature

• Timer starts when the injection is initiated. This feature is especially useful to time the start of a delayed phase scan.

Scan Delay

• Enables preset or injection delay to maximize synchronization with the scanner.

Multi-phase Option

• Up to six phases and Hold & Pause options permit full capability for today's complicated scan protocols.

Protocol Lock

• Current protocol can be locked at the display to enable remote arming and injections from the patient's side to enhance safety.

medrad











Medrad is dedicated to continually improving the quality of all our products and services such that our customers' satisfaction, loyalty, and respect are unsurpassed. It is our policy to clearly understand and agree upon the valid requirements of the work we perform for our customers, both internal and external, and to pursue 100% conformance to those requirements... on

time, every time. Medrad reserves the right to modify the specifications and features described herein, or discontinue manufacture of the product described at any time without prior notice or obligation. Please contact your authorized Medrad representative for the most current information.

Medrad. Inc. is ISO 9001/ISO 13485 Certified.

Medrad, Quality for Life, Q and design are federally registered trademarks, Stellant is a trademark, and Medrad is a service mark of Medrad, Inc. U.S.A.





Injection Specifications

Miscellaneous Specifications	
Maximum Number of Protocols	32
Maximum Number of Phases	6
Syringes (Volume capacity)	200 ml sterile disposable syringe
Hold	Maximum HOLD time is 20 minutes
Pause	1-900 seconds (15 minutes) in 1 second increments
Scan Delay	0-300 seconds (5 minutes) in 1 second increments
Programmable Pressure Limit psi/ kPa	200 ml syringe: 325 psi, 2241 kPa
Volume (range & increments)	1 ml to syringe capacity in 1 ml increments
Flow Rate (range & increments)	0.1 to 10 ml/sec in 0.1 ml increments

Electrical Requirements (VAC/Hz)	100-240 VAC, 50/60 Hz
Syringe Heater Range	95 degrees F +/-9 degrees 35 degrees C +/-5 degrees
Head Mounting Options	Overhead Counterpoise System or Floor Pedestal (Floor Pedestal includes an Integral IV Pole)

Mechanical				
Specifications	Height	Width	Depth	Weight
Single Injector Head	6.1" (15,5 cm)	8.6" (21,8 cm)	14.5" (36,8 cm)	12 lbs (5,4 kg) (without syringe)
Display Control Unit with stand	13.5" (34,3 cm)	12.5" (31,9 cm)	8.9" (22,5 cm) Screen only = 3.4" (8.6 cm)	8.2 lbs (3,7 kg)
Base Unit	11.5" (29,2 cm)	11.0" (27,9 cm)	8.8" (22,2 cm)	13.6 lbs (6,2 kg)

Ordering Information

Systems listed include all standard cabling

Catalog Number: SCT 121 Stellant Sx Injection System with Pedestal Head Mount

Catalog Number: SCT 122 Stellant Sx Injection System with Counterpoise System



Quality, Medrad Service and Medrad Support.

© 2003 Medrad, Inc. All Rights Reserved.

STL EN SX 202021

Medrad, Inc. Corporate Office One Medrad Drive Indianola, PA 15051-0780 USA 412-767-2400

www.medrad.com For more information: info@medrad.com

Customer Service/Orders 1-800-MEDRAD-1 (1-800-633-7231)

Customer Service FAX (412) 767-4120

International FAX (412) 767-4128