


External Dimensions (Unit:mm)



Founded in 1875, Shimadzu Corporation, a leader in the development of advanced technologies, has a distinguished history of innovation built on the foundation of contributing to society through science and technology. We maintain a global network of sales, service, technical support and applications centers on six continents, and have established long-term relationships with a host of highly trained distributors located in over 100 countries. For information about Shimadzu, and to contact your local office, please visit our Web site at www.shimadzu.com



SHIMADZU CORPORATION. International Marketing Division 3. Kanda-Nishikicho 1-chome, Chiyoda-ku, Tokyo 101-8448, Japan Phone: 81(3)3219-5641 Fax. 81(3)3219-5710 URL http://www.shimadzu.com




Shimadzu Corporation Medical Systems Group has been certified by TÜV Rheinland as a manufacturer of medical equipment and systems in compliance with ISO9001:2000 Quality Management Systems and ISO13485:2003 Medical Equipment Quality Management Systems.


Remarks;

- *Every value in this catalogue is a standard value, and it may vary a little from the actual at each site.
- a little from the actual at each site.

 *The appearances and specifications are subject to change for reasons of improvement without notice.

 *Certain configurations may not be available pending regulatory clearance. Contact your Shimadzu representative for information. on specific configurations.
- *Before operating this system, you should first throughly review the Instruction Manual.







VB9 Slender



Direct Conversion Provides Outstanding Image Sharpness

Fast, Easy, and Precise Positioning

Leading-Edge 3D Application Provides a Valuable IVR Tool

Carefully Designed Technology to Help Reduce Exposure Levels to Patients and Medical Staff

System Design Improves Workflow Efficiency

State-of-the-Art Digital Technology Ensures Faster Response

Shimadzu, a pioneer in angiography systems using direct-conversion FPDs, has developed a bi-plane system that produces even higher quality images and ensures interventional procedures are safe and smooth.

The ultra-high-definition images provided by safire's direct conversion technology are outstanding for rendering fine blood vessels in the head, as well as showing clearly all guide wires, microcatheters, coils, and other sensitive devices.

And installation of our new image processor SUREengine has made it possible to improve image quality even further.

Perform advanced-level treatments (interventions) using slender C-arms that, together with this system's ultra-compact design, ensure easy positioning during bi-plane operation, together with a multiprocessor-equipped digital system for fast response.



Flat Panel Detector Direct Conversion Provides Outstanding Image Sharpness





clearly, the design of this system's directconversion FPD uses extremely small pixels

of just 150 um.

BRANSIST safire's clear rendering of microvessels and sensitive devices used for interventional procedures provides a powerful tool for examinations and treatments.







Ultra-High-Definition Real-Time DSA

A direct-conversion FPD makes ultra-high-definition real-time DSA possible and can be applied not only to intervention procedures in head regions, but the entire body as well, including the chest and extremities.






Compared to images produced by other angiography systems, BRANSIST safire's direct-conversion FPD delivers the highest quality.


SUREengine, our real-time image processor designed specifically to enhance image quality, performs real-time contrast enhancement and noise reduction on the respective frequency bands of fluoroscopic and radiographic images, ensuring improved visibility of intervention devices and blood vessels.




Equipped with a Next-Generation X-Ray Detector to Create the Ideal Imaging Chain

Direct-conversion detection, which converts X-rays directly to electrical signals, eliminates the blurring that results from the scattered light generated in conventional imaging processes, realizing "pure digital imaging". It also features the smallest pixel size in its class – $150 \mu m$ – to provide fluoroscopy/radiography images with exceptional visibility.









Intervention Support Fast, Easy, and Precise Positioning

Compared to single-plane systems, bi-plane system C-arm positioning and entry of settings are more complicated.

However, BRANSIST safire is equipped with proprietary features such as CyberChase, CyberGrip, DirectMemory, and TriplePivot, which ensure fast C-arm positioning.



Versatile, One-Hand Control of the C-Arm via CyberGrip

The CyberGrip one-hand controller allows the X-ray technologists to move the high-speed C-arm freely.

This controller is equipped with a grip-type handle that provides intuitive operation so that the C-arm, table, and FPD can all be moved without changing your grip.

This frees X-ray technologists to concentrate on the actual procedures, without having to worry about the position of their hand.




DirectMemo

New Concept Auto Positioning [DirectMemory]


The new DirectMemory unit intuitively calls up C-arm clinical angles for a patient and allows rapid repositioning.

DirectMemory provides powerful support for smooth, stress-free diagnostic catheterization and interventional procedures.




High-Speed C-Arms Ensure Quick Examinations

The C-arms are capable of single-plane positioning at 25°/s (frontal arm) and bi-plane positioning at 15°/s. You can also quickly switch between single-plane motion and bi-plane motion.



CyberChase Automatically Tracks the Area of Interest

The imaging chain for the lateral C-arm moves vertically up to 17.5 cm, allowing you to quickly adjust the bi-plane position to a new area of interest, without changing table height. Even if the C-arm angle is changed, CyberChase automatically tracks the area of interest and allows even complicated C-arm positioning to be performed quickly and accurately.





New Application

Leading-Edge 3D Application Provides a Valuable IVR Tool

BRANSIST safire includes state-of-the-art 3D imaging software, effective for complicated interventions that require super-selective catheterization and wide-ranging therapy. Thanks to this, you can perform high-level medical procedures reliably, safely, and efficiently.

3D Angiography Safire 3D-A (optional)

The combination of ultrahigh-definition rotational DA and DSA images obtained via direct-conversion FPDs and a state-of-the-art 3D reconstruction workstation (NaviDAS+ 3D-Recon) equipped with a wide range of functions ensure that this option produces superior diagnostic images. Additionally, our Max Scan function quickly rotates the C-arm at 60°/second, shortening the time required for large-range radiography. This in turn reduces injection time and usage volume of contrast medium, helping to lower patient stress.






New Standard-Equipment Motion-Tolerant RSM-DSA System


Based on Shimadzu high-speed digital image processing technology, this new RSM-DSA system provides extremely motion-tolerant DSA without requiring a mask run, by subtracting digital frequency processed mask images from live images.

This makes Shimadzu's RSM-DSA especially effective when performing contrast examinations on patients who have difficulty holding their breath, and for 3D contrast examinations that combine bolus chase examinations of lower limbs with precession, pendular, or other C-arm movements.









Carefully Designed Technology to Help Reduce Exposure Levels to Patients and Medical Staff

With treatments involving intervention procedures on the increase, exposure levels to physicians, X-ray technologists and patients cannot be ignored.

The BRANSIST safire system does not rely merely on its high-performance direct-conversion FPD to determine a trade-off between low exposure and high image quality, but rather accomplishes a balance using the entire system.


Built-In MBH Filter Excludes Soft X-Rays

BRANSIST safire has MBH (Multi Beam Hardenina) filters built into the front end of the X-ray tube.

Automatically select the optimum filter for examination or operation and to exclude soft X-rays that do not contribute to the image.

Wave-Tail Elimination via Grid Control

The X-ray generator uses grid control to enable the ideal form of wave-tail elimination. In low-dose pulsed fluoroscopy mode, accurately eliminating X-rays corresponding to wave tails helps reduce exposure even further.





Reducing the Operator's X-Ray Exposure **LL**↔**RL** Positioning

One touch of a button is all that is needed to position the lateral arm's X-ray tube on the opposite side of the physician (R-L position) so as to minimize the scattered radiation.

"It meets ICRP* Pub.85" *ICRP:International Commission on Radiological Protection.

Overall **Reduction in Exposure**









Workflo

System Design Improves Workflow Efficiency

Equipped with state-of-the-art digital technology and a parallel microprocessor, BRANSIST safire systems can simultaneously process fluoroscopy and radiography images.

Together with a large-capacity tube and larger image storage capacity, this facilitates a more efficient workflow.

Equipped with High-Capacity Liquid-Bearing X-Ray Tube

A 3.0-MHU, high-capacity, liquid-bearing X-ray tube (LX-3081) is a standard feature, provides a quiet examination room and eliminates the need for cooling time between repeated examinations.

This allows efficient operation at facilities performing large numbers of examinations.

Multiprocessor for Parallel Processing

The system includes a multiprocessor capable of processing image data, such as archiving images, sending images via a network, or processing or measuring images, even while collecting radiographic or fluoroscopic images.


High-Speed Image Acquisition

BRANSIST safire can acquire images during bi-plane operation at up to 30 fps using alternating exposures.



Seamless Interface with DICOM Network

BRANSIST safire provides flexible handling of image networking processes (DICOM storage, storage commitment, Query/Retrieve services) and HIS management (optional DICOM MWM service), making it possible to interface seamlessly with the facility network.





Fast Respons

State-of-the-Art Digital Technology Ensures Faster Response

The speed at which X-ray technologists can display necessary or desired images is an extremely important factor in performing neurointervention procedures that are free of stress.

Using state-of-the-art digital technology and multiprocessors, BRANSIST safire can instantly select and display both still and moving reference images.

Also, multiprocessor technology provides the benefit of improved image data security through the use of parallel image processing and data mirroring.

Fast Bedside Imaging

The joystick type IVR NEO bedside image controller includes only the features that are essential for intervention procedures. Thanks to outstanding responsiveness in reference image selection, multi-image display, fluoroscopy/radiography mode switching, map mode selection, fluoroscopy log, and high speed replay, X-ray technologists are free to perform procedures in any manner they desire.



Fluoro Back Record

This feature keeps a record of scenes up to 1023 frames back and allows retroactively recording any of these scenes onto the hard drive.

This is especially useful for recording the important fluoroscopy scenes or MAP images.

These images can also be stored as DICOM images on external media or transferred to other locations via a network.



Data Mirroring

By installing two hard drives for storing data, images are instantaneously mirrored to a second computer, which allows not only processing in parallel, but also provides a backup for peace of mind.



Dual Road Mapping

Standard features of this system include two types of road mapping functions. One type is a superimposed mapping of a fluoroscopic image overlaid on a DSA image.

The other type is a subtraction mapping where a peak-hold image is flashed during fluoroscopy processing and then subtracted from the live fluoroscopic image. This allows easy selection of the optimal mapping mode for a specific technique or objective.