



Founded in 1875, Shimadzu Corporation, a leader in the development of advanced technologies, has a distinguished history of innovation built on the foundation of contributing to society through science and technology. We maintain a global network of sales, service, technical support and applications centers on six continents, and have established long-term relationships with a host of highly trained distributors located in over 100 countries. For information about Shimadzu, and to contact your local office, please visit our Web site at www.shimadzu.com



SHIMADZU CORPORATION. International Marketing Division 3. Kanda-Nishikicho 1-chome, Chiyoda-ku, Tokyo 101-8448, Japan Phone: 81(3)3219-5641 Fax. 81(3)3219-5710 URL http://www.shimadzu.com



Remarks;

- *Every value in this catalogue is a standard value, and it may vary a little from the actual at each site.
- *The appearances and specifications are subject to change for reasons of improvement without notice.
- *Certain configurations may not be available pending regulatory clearance. Contact your Shimadzu representative for information on specific configurations.
- on specific configurations. %Before operating this system, you should first throughly review the Instruction Manual.




C506-E051C





High Quality and Large Field of View FPD Imaging





12 million Add

17-inch FPD Covers a Large View Field

The 17 x 17-inch maximum view field applies the extremely high definition and density resolution of the direct-conversion FPD across an area equivalent to 14 x 17-inch film. The large distortion-free view field means only minimal positioning movements are required. This ensures trouble-free handling of diverse examination regions, including orthopedic procedures and abdominal imaging of the stomach and large intestine.


High-Definition, Moving Images

Shimadzu's direct-conversion FPD is able to capture ultra-fine targets in real time. Images that were difficult to display on previous systems can now be reproduced with great accuracy


View field size	43 x 43 cm
X-ray conversion method	Direct-conversion a-Se
Pixel pitch	150 µ m
Resolution (max.)	3.3 lp/mm
Dynamic range	14 bit
Frame rate	30 fps

Equipped With a Next-Generation X-Ray Detector to Create the Ideal Imaging Chain

The direct-conversion method, where X-ray input signals are directly converted to electric signals, completely eliminates the intermediate noise that occurs with conventional imaging processes. This realizes true, pure digital imaging based on the linear conversion of X-ray signals to electric signals.




Wide-View 17-inch FPD



Dynamic Image Observations Over a Large View Field

During contrast medium examinations, the FPD's wide view field offers distortion-free dynamic images without moving the imaging system or tabletop.

This provides greater clarity when observing the flow of barium from the esophagus to the stomach



Esophagus Serial Images Using Wide View Field FPD

Distortion-Free, High-Definition Images Over a Wide View Field

Unlike an image intensifier, the flat-panel detector offers a planar and rectangular detector face. The net effect of this technology is uniform resolution and contrast from the center right to the edges. The FPD is unaffected by magnetic fields, ensuring distortion-free images. These high-definition images permit highly accurate measurements even if the target is at the periphery.





Image Intensifier (Comparison using Shimadzu products)

Wide View Field and Long Stroke **Cover a Large Examination Area**

cover an extensive imaging area.






For Bariatric Imaging

The table's heavy-duty design supports a patient load of up to 318 kg (700 lb) in the horizontal position. This is the best in its class

The Wide Field of View **Tomography Provides** Accurate Images of Area of Interest

Tomography is possible with the patient standing or in a reverseinclined position. Tomography in a loaded state gives you a lots of clinical benefits which are not visualized by the conventional horizontal tomography.

The FPD's wide 17-inch view field combines with the large longitudinal stroke of the Sonialvision system to



Lift Mechanism

Shimadzu's proprietary table-lift mechanism has been designed for the comfort of both the patient and operator

The tabletop lowers to just 47 cm to allow radiography procedures with minimum stress on the patient. It is also easy to adjust the table to a convenient height for patient transfer or surgical procedures.

.1.5 m

Radiography at 1.5 m Accommodates Patients in Wheelchairs

The system is able to maintain an adequate distance to accommodate radiography examinations on patients in wheelchairs.

Easy Examination of Knee Joints and Leg Veins

A full stroke of the imaging system is possible at all table tilt angles. This ensures safety during examinations such as knee-joint kymography and leg venography.

Upgraded Digital Unit Fully Utilizes the Direct-Conversion FPD

Only Shimadzu's direct-conversion FPD is capable of handling 2880 x 2880-matrix ultra-high-definition images. A cutting-edge digital system directly processes and stores these high-precision images at the original capture quality.

Simple Windows[®]-Based Operability

This new digital radiography unit is designed around the highly reliable Windows[®] operating system, which makes the rapid processing of high-definition images possible using simple operations. The graphical user interface and mouse provide an intuitive operation environment that is easy to learn and easy to use.

Large Capacity Digital Image Storage

Accumulated digital images are recorded in real-time to an internal large-capacity high-speed hard disk.

Images can also be saved externally to DVD-R/CD-R media, allowing storage of up to 2000 DICOM format frames on a single disk.




2880 x 2880 Matrix **Original Quality Image Acquisition**

The 2880 x 2880-matrix ultra-high-definition digital images generated by Shimadzu's direct-conversion FPD permit the observation of extremely fine structures. When stored, each radiographic image maintains the original quality of the matrix at the time of capture. The images retain their outstanding quality even when enlarged on a viewer.



High-Speed Imaging for Perfect Radiography Timing


The direct-conversion FPD is able to capture images at high speed without ghosting or burn out. This makes it extremely valuable for gastrointestinal and other angiographic examinations.



High-Speed High-Definition Radiography Up to **15** frames per second

Programmable Positioning Switches Color LCD Panel

Examination start positions, such as table height, and post-examination return positions can be programmed to accommodate the operating environment. This customization permits faster turnaround on procedures.



2

201327





Upper Gastrointestinal Mode Radiography programs can be set for spot and serial radiography.



X-ray setting screens. The design of these easy-to-recognize color-coded screens



Multi-Frame Imaging Using Similar **Operations With Film Radiography**

Digital multi-frame imaging can be performed using similar operations to multiframe imaging with conventional film radiography.

The multi-frame programs available include horizontal split, vertical split and four-way split.

Technology That Expands Your Application Possibilities

A full and clear understanding of customer medical imaging expectations and needs enables Shimadzu to continuously provide a range of effective solutions. These solutions are the result of creative, flexible thinking and provide even greater clinical benefits when coupled with the long stroke and high degree of freedom of our Sonialvision table, and the ultra-wide field-of-view and high-resolution images of our 17-inch FPD.

17-Inch Wide Field-of-View DSA

Motion Tolerant RSM-DSA OPTION

Breath-Holding Not Required

Freedom in Framing

medium

gases or patient body movements.

No Need For Mask Image Reduces Exposure

High-speed data acquisition

MAX 7.5 fps

No mask run is required, allowing use of far less contrast medium.

images.

Shimadzu's proprietary Real-Time Smoothed Mask Digital Subtraction Angiography (RSM-

DSA) is a revolutionary new DSA application that eliminates the need for acquiring mask

Provides higher diagnostic-value DSA images with no artifacts from digestive tract

No artifacts from body movements or breathing eliminates need to immobilize the patient and allows movement of imaging location together with flow of contrast

Shimadzu's direct-conversion FPD is able to capture ultra-fine targets in real time. Images that were difficult to display on previous systems can now be reproduced with great accuracy.



RSM Filter

An advancement in RSM-DSA technology, the new high-speed digital image processing filter provides real time viewing of fluoroscopy and radiography images. This filter's production of a more uniform background suppresses halation and provides images with higher visibility.

Tomosynthesis (Digital Multislice Tomography) OPTION

Tomosynthesis uses digital data from a single tomographic scan to create an image of slice desired

Expands Examination Range

Tomosynthesis allows recording of images at any angle required for diagnosis, including a variety of table angles or with the patient standing to apply gravity.

Reduces Examination Time and Eliminates Need For Additional Imaging

Since any slice can be reconstructed from a single tomographic scan, total recording time is much less than conventional digital tomography.

Also, fewer imaging errors help to shorten patient exposure and immobilization times, reducing patient stress.

Fewer Metal Artifacts

- Useful for Orthopedic Examinations

Slice images can be viewed without the influence of metal artifacts that is often seen in CT images. This is useful in follow-up examinations after orthopedic surgeries where metal objects have been implanted.



TOMOSYNTHESIS

RSM-DSA

Slot Radiography (Longitudinal Exposure) OPTION

Parallel movement of the imaging chain and FPD unit with a slit-collimated X-ray exposure is used to produce a longitudinal radiographic image. Wide-range longitudinal images are then

acquired by stitching of individual slit images Shimadzu's unique 1.5 m imaging chain

extension enables smooth radiographic image processing of the spine and lower extremities in both the standing and supine position.

scatter.

made with the system.

DO BHIMADZU

Dual Energy Subtraction (D. E.) OPTION

Energy subtraction between high- and low-voltage images enables separate display of soft tissue and bone images. This allows diagnosis of nodular density hidden behind the ribs in soft tissue images, and of calcification in bone images. Diagnostic D.E. is possible using the serial radiography function at a maximum of 15 fps. This enables diagnostic imaging of the soft tissue of the lungs while reducing visible respiratory movement.

DUAL ENERGY SUBTRACTION

High-Quality Longitudinal Images

In addition to the extremely high-definition image acquisition of the direct-conversion FPD, highquality longitudinal images are processed using slit X-ray exposures, while minimizing X-ray

Minimal-Distortion Images Ensure **Highly Accurate Measurement Functions**

This system incorporates measurement functions for distance, Cobb angle and others.

In addition to the distortion-free images produced by the FPD, the effect of oblique X-rays is minimized by the slit-collimated X-ray exposures. This allows highly accurate measurements to be



SLOT RADIOGRAPHY

Tomosynthesis Dual Energy Subtraction OPTION

Combination of Dual Energy Subtraction with Shimadzu's Tomosynthesis function enables the acquisition of 3D information of soft tissue and bone.

This allows 3D observation of bronchial tubes and nodular density in the lung and of articular cavities by subtracting the bone image.

TOMOSYNTHESIS

Comprehensive Dose Management for Patient and Operator

The increasing use of IVR treatment requires a system that can offer both high image quality and low exposure dose. The Sonialvision Safire II with direct-conversion FPD meets this requirement by incorporating an exposure reduction mechanism matched to the high image performance of the FPD. These two technologies work in tandem to achieve an overall balance between high image quality and low exposure.

Low-Dose Pulsed Fluoroscopy


Four low-dose pulsed fluoroscopy modes are available. They include the wave-tail cut-off mechanism, which obtains a more ideal pulsed waveform to achieve further exposure reductions.

Pulsed Fluoroscopy further.



Iris Collimator Promotes Image Quality and Exposure Reduction

A precision iris collimator is used to efficiently cut out unnecessary areas. Collimating the required view field for IVR or repositioning effectively restricts image deterioration due to scattered X-rays and also reduces the dose by cutting out unnecessary X-rays.




Conventional rectangular collimator


Ideal Wave-Tail Elimination Grid Control

The X-ray generator uses grid control to enable the ideal form of wave-tail elimination. In low-dose pulsed fluoroscopy mode, accurately eliminating X-rays corresponding reduce exposure during fluoroscopic to wave tails helps reduce exposure even

Digital Recording of Fluoroscopic Images During Examinations

Fluoroscopic images can be recorded in the internal memory at up to 30 fps during examinations. These images can be viewed using instantaneous or cyclic playback to diagnosis.








BH Filter Eliminates Superfluous X-Rays


The system incorporates a BH (Beam Hardening) filter to eliminate soft X-rays. This BH screening cuts out X-rays that do not contribute to the image



Safe and Comfortable Examination Environment


Contact Safety Switch

The imaging system is equipped with a safety switch that automatically halts movement in the unlikely event the equipment contacts an object during a procedure. This important feature ensures the safety of both the operator and patient.




Overhang Detector Switch

The imaging system also features a second safety switch that automatically stops movement if a finger or other body part overhangs the table during an examination. This switch offers an extra safeguard to the patient.



Total Network Solutions for the High-Speed IT Age



Side Station OPTION

The side station enables the use of DICOM-compatible viewing Remote service support is readily available by connecting the system to allow DR image observation, processing, printout and digital image-processing unit to a Shimadzu maintenance facility server transfer to be performed during examinations. by standard telephone line. It also enables the use of the tomosynthesis feature that Shimadzu offers rapid system diagnosis and remedy in the unlikely

reconstructs images of any arbitrary plane. For Sonialvision Safire Il use only.

Configuration	-		· · · · · · · · · · · · · · · · · · ·		1
Remote-operated R/F table	ZS-100I		Diagnostic X-ray high-voltage unit X-ray tube unit		UD150B-40 (80 kW)
Digital radiography unit	DAR-8000	of (including FPD)			0.7/1.2JG326D-265AT (750 kl
Local-controrol console	1	Side station	II Special Accessories	Tomosy	
Local-controrol console Shoulder rest	1	Side station Endoscope holder		Tomosy RSM-D	
				RSM-D	
Shoulder rest	1	Endoscope holder		RSM-D Slot Rad	SA
Shoulder rest Hand grip (two types)	1 1 set	Endoscope holder Knee clutches		RSM-D Slot Rad Dual En	SA diography
Shoulder rest Hand grip (two types) Table mattress	1 1 set 1	Endoscope holder Knee clutches Drain bag		RSM-D Slot Rad Dual En	SA diography ergy Subtraction storage
Shoulder rest Hand grip (two types) Table mattress Compression cone (flat and extruded)	1 1 set 1 1 set	Endoscope holder Knee clutches Drain bag Cystographic chair		RSM-D Slot Rad Dual En DICOM	SA diography ergy Subtraction storage MWM

Remote Maintenance

event of a problem.