



CardioMD™

compact, modern, and dedicated cardiac camera

The CardioMD is a fixed 90 degree dedicated nuclear cardiology camera specifically designed to meet the unique requirements of an office-based practice. With excellent image quality, high throughput, and a small footprint, the CardioMD offers complete cardiac imaging capabilities. This innovative system provides the best solution for optimal departmental workflow with easy to use features and simple patient setup. The CardioMD is packaged with state-of-the-art cardiology software.



FEATURES

- Compact camera
- Open gantry design
- Unique table design
- Windows-based workstation
- Acquisition display capability

BENEFITS


- Optimized for the office environment, fitting into a room as small as 8' x 10'
- Fast patient setup with complete access during exam
- Unobstructed cardiac imaging with 4" table cut-offs
- Enhanced workflow with effortless data access
- Ability to review and QC the "just" acquired images on the acquisition station in a cine or frame mode display




Acquisition Workflow




Step 1: Setup Patient



Step 2: Acquire Study



Step 3: Display Images



The CardioMD offers the appeal of a lightweight, comprehensive cardiac camera with the convenience of a two-day installation period. CardioMD's acquisition station guides technologists through the daily workflow with simultaneous display of both detectors, spectrum p-scope, local data storage, and immediate display capability. The acquired images are processed using state-of-the-art software including AutoSPECT™Plus, AutoQUANT™, CardiaQ™, and Emory Toolbox™.

PHILIPS MEDICAL SYSTEMS

540 Alder Drive
Milpitas, CA 95035
800-538-8351
www.medical.philips.com

© 2003 Philips Medical Systems Nederland B.V. ADAC, the ADAC logo, GEMINI, ALLEGRO, and PIXELAR are trademarks or registered trademarks of ADAC Laboratories, Inc. All other trademarks or registered trademarks are property of their respective holders. All rights are reserved. Reproduction in whole or in part is prohibited without the prior written consent of the copyright holder.

Philips Medical Systems Nederland B.V. reserves the right to make changes in specifications or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

4535 983 03198/882 2003-05

